

Tekniske data

Fluke 700HTPK Hydraulic Test Pressure Kit



Nøglefunktioner

- Easy to use, rugged and reliable construction, capable of generating hydraulic pressure up to 10,000 psi/690 bar
- Combine with one of the Fluke 700G Precision Pressure Gauges for a complete pressure testing and calibration solution
- In addition, combine with the Fluke 700G/TRACK Software to upload over 8,000 pressure measurements logged in the field to a PC

Produktoversigt: Fluke 700HTPK Hydraulic Test Pressure Kit

Generates pressure up to 10,000 psi, 690 BAR. Includes the 700HTP-2 hydraulic hand pump, test hoses and adapters, and protective hard case.

Modeller



Fluke-700HTPK

Hydraulisk testtryksæt med tryktestmåling med god præcision, hydraulisk håndpumpe, 1 slangesæt og adaptere, indbefatter hård kuffert

Includes:

- 1 hose kit and adapters
 - Hard case
-

Fluke. *Keeping your world up and running.®*

Fluke Danmark A/S
c/o Radiometer Medical ApS
Åkandevej 21
2700 Brønshøj
Danmark
Tlf.: 70 23 58 53
E-mail: cs.dk@fluke.com
www.fluke.dk

©2022 Fluke Corporation. Alle rettigheder
forbeholdes.
Oplysningerne kan ændres uden forudgående
varsel.
01/2022

**Ændringer i dette dokument er ikke tilladt uden
skriftlig tilladelse fra Fluke Corporation.**