

Tekniske data

Fluke 700TC1 Thermocouple Plug Kit



Produktoversigt: Fluke 700TC1 Thermocouple Plug Kit

Et sæt med ti mini-stikkonnetorer:

Modeller



Fluke 700TC1

Fluke 700TC1 Thermocouple Plug Kit
Thermocouple Plug Kit

A kit of 10 mini-plug connectors

- Type J (black), one
- Type K (yellow), one
- Type T (blue), one
- Type E (purple), one
- Type R/S (green), one
- Type B or CU (white), one
- Type L (J-DIN) (blue), one
- Type U (T-DIN) (brown), one
- Type C (red), one
- Type N (orange), one
- One year warranty

Fluke. *Keeping your world up and running.®*

Fluke Danmark A/S
c/o Radiometer Medical ApS
Åkandevej 21
2700 Brønshøj
Danmark
Tlf.: 70 23 58 53
E-mail: cs.dk@fluke.com
www.fluke.dk

©2021 Fluke Corporation. Alle rettigheder
forbeholdes.
Oplysningerne kan ændres uden forudgående
varsel.
12/2021

**Ændringer i dette dokument er ikke tilladt uden
skriftlig tilladelse fra Fluke Corporation.**