


Tekniske data

## Fluke ES-162P4-2 4-Pole Stake Kit



### Produktoversigt: Fluke ES-162P4-2 4-Pole Stake Kit

Needed to perform 3- and 4-Pole Fall-of-Potential testing (Stake testing). Used with both the Fluke 1623-2 and Fluke 1625-2. Already included in the Fluke 1625-2 Kit.

## Modeller



### Fluke ES-162P4-2

- (4) Stakes
- (1) 50 m cable reel of red wire
- (1) 25 m cable reel of blue wire
- (1) 25 m cable reel of green wire

**Fluke.** *Keeping your world up and running.®*

**Fluke Danmark A/S**  
c/o Radiometer Medical ApS  
Åkandevej 21  
2700 Brønshøj  
Danmark  
Tlf.: 70 23 58 53  
E-mail: [cs.dk@fluke.com](mailto:cs.dk@fluke.com)  
[www.fluke.dk](http://www.fluke.dk)

©2021 Fluke Corporation. Alle rettigheder  
forbeholdes.  
Oplysningerne kan ændres uden forudgående  
varsel.  
12/2021

**Ændringer i dette dokument er ikke tilladt uden  
skriftlig tilladelse fra Fluke Corporation.**