

Tekniske data

VPS200-G Voltage Probe Set, Gray 200 MHz (one gray)



Produktoversigt: VPS200-G Voltage Probe Set, Gray 200 MHz (one gray)

Double insulated voltage probe for ScopeMeter190 series.

Specifikationer: VPS200-G Voltage Probe Set, Gray 200 MHz (one gray)

Specifications	
Attenuation	10x
Bandwidth	200 MHz
Max. input voltage	1000 V
Max. voltage to ground	600 V CAT III 1000 V CAT II per IEC / EN61010-031
Input resistance	10 Mohm
Input capacitance	14 pF
Compensation range	10 to 22 pF



Cable length	1.2m
Warranty	1 Year
Specifications apply when the probe is connected to Fluke 190 Series Scopemeter.	

Modeller



Fluke VPS200 G

Voltage Probe Set, Gray 200 MHz (one gray)

Each probe set includes:

- 600 V CAT III / 1000 V CAT II rated:
 - 10:1 Voltage Probe
 - Hook Clip for Probe Tip, 4-mm Test Probe for Probe
 - Ground Lead with Hook Clip
 - Ground Lead with Mini Alligator
 - The red probe set also includes a 30 V CAT II rated Ground Spring for Probe Tip
-

Fluke. *Keeping your world up and running.®*

Fluke Danmark A/S
c/o Radiometer Medical ApS
Åkandevej 21
2700 Brønshøj
Danmark
Tlf.: 70 23 58 53
E-mail: cs.dk@fluke.com
www.fluke.dk

©2021 Fluke Corporation. Alle rettigheder
forbeholdes.
Oplysningerne kan ændres uden forudgående
varsel.
11/2021

**Ændringer i dette dokument er ikke tilladt uden
skriftlig tilladelse fra Fluke Corporation.**