

Tekniske data

700PMP Trykpumpe



Produktoversigt: 700PMP Trykpumpe

Description: The Fluke 700PMP is a hand-operated pressure pump to provide pressures up to 150 psi/1000 kPa. Output fitting is 1/8 FNPT.

Application: Linear stroke of 4 cm. Multi-turn vernier for fine adjustment of pressure.

Specifikationer: 700PMP Trykpumpe

Mechanical and General Specifications

Size	265 x 35 mm
Weight	340 g
Warranty	One year

Modeller



Fluke 700PMP

PMP Pressure Pump

Instruction sheet

Fluke. *Keeping your world up and running.®*

Fluke Danmark A/S
c/o Radiometer Medical ApS
Åkandevej 21
2700 Brønshøj
Danmark
Tlf.: 70 23 58 53
E-mail: cs.dk@fluke.com
www.fluke.dk

©2021 Fluke Corporation. Alle rettigheder
forbeholdes.
Oplysningerne kan ændres uden forudgående
varsel.
12/2021

**Ændringer i dette dokument er ikke tilladt uden
skriftlig tilladelse fra Fluke Corporation.**