

Tekniske data

Fluke CO-205 Aspirator Kit



Nøglefunktioner

- Stainless steel sampling tube
- Industrial-grade hand operated aspirator to draw flue sample
- Easy to replace particulate filter
- Specially designed nose cap for connection to the Fluke CO-210/220
- 4' (1.2 m) hose length

Produktoversigt: Fluke CO-205 Aspirator Kit

Sample flue gas with your Fluke CO-210 or CO-220 using the CO-205 Aspirator Kit

Draw flue gas samples up to 700°F (370°C) reliably for measurement with the Fluke CO-210 Carbon Monoxide probe or Fluke CO-220 Carbon Monoxide meter.

The CO-205 aspirator accessory kit contains all the components necessary to provide a clean sample for the Fluke family of gas measuring devices.

Modeller



Fluke CO-205

Fluke CO-205 Aspirator Kit

Fluke. *Keeping your world up and running.®*

Fluke Danmark A/S
c/o Radiometer Medical ApS
Åkandevej 21
2700 Brønshøj
Danmark
Tlf.: 70 23 58 53
E-mail: cs.dk@fluke.com
www.fluke.dk

©2021 Fluke Corporation. Alle rettigheder
forbeholdes.
Oplysningerne kan ændres uden forudgående
varsel.
12/2021

**Ændringer i dette dokument er ikke tilladt uden
skriftlig tilladelse fra Fluke Corporation.**