

Tekniske data

P5544 Two Gauge Stand



Produktoversigt: P5544 Two Gauge Stand

This adapter mounts directly to the test port of the calibrator and allows for the calibration of two instruments at the same time or the connection of a reference test instrument. Maximum working pressure is 10,000 psi (700 bar).

Specifikationer: P5544 Two Gauge Stand

Specifications	
Body Material	Brass and Stainless steel
O-Rings	Nitrile
Test Port Adapters	1 set, 1/8, 1/4, 3/8 and 1/2 NPT or BSP and adapters supplied with calibrator
Calibrator Test Connection	1/2 BSP male (mating adapter supplied)

Modeller



P5544

Two Gauge Stand

Two Gauge Stand

Fluke. *Keeping your world up and running.®*

Fluke Danmark A/S
c/o Radiometer Medical ApS
Åkandevej 21
2700 Brønshøj
Danmark
Tlf.: 70 23 58 53
E-mail: cs.dk@fluke.com
www.fluke.dk

©2021 Fluke Corporation. Alle rettigheder
forbeholdes.
Oplysningerne kan ændres uden forudgående
varsel.
11/2021

**Ændringer i dette dokument er ikke tilladt uden
skriftlig tilladelse fra Fluke Corporation.**