

Elektrotester-Kit, T5-1000 mit Holster und 1AC II



Wichtigste Merkmale

- T5-1000 electrical tester automatically measures volts AC and volts DC with precise digital resolution
- Includes easy, accurate OpenJaw™ current measurement
- Displays resistance to 1000 Ω
- Includes a continuity beeper for quick go/no-go tests
- All inputs are protected to EN61010-1 CAT III 1000 V CAT IV 600 V. UL, CSA
- 1AC II offers Voltbeat™ technology and continuous self-test so you always know it's working
- Upon detection of voltage, the VoltAlert™ tip glows and beeper sounds
- CAT IV 1000 V for added protection

Produktübersicht: Elektrotester-Kit, T5-1000 mit Holster und 1AC II

Spannungs-, Durchgangs- und Stromtesterkit mit Holster. Der 1AC ist kostenlos.

Modelle



Fluke T5-H5-1AC II Kit

This kit includes:

- Fluke T5-1000
 - Fluke 1AC II
 - H5
-

Fluke. Damit Ihre Welt intakt bleibt.

Fluke (Switzerland) GmbH

Industrial Division
Hardstrasse 20
CH-8303 Bassersdorf
Telefon: +41 (0) 44 580 7504
E-Mail: info@ch.fluke.nl
www.fluke.ch

©2021 Fluke Corporation. Alle Rechte vorbehalten.
Anderungen vorbehalten
11/2021

**Dieses Dokument darf nicht ohne die schriftliche
Genehmigung der Fluke Corporation geändert
werden.**

Technischer Beratung:

Beratung zu Produkteigenschaften, Spezifikationen,
Messgeräte und Anwendungsfragen
Tel.: +49 (0) 7684 8 00 95 45
E-Mail: techsupport.dach@fluke.com