

TECHNICAL DATA

## Fluke BP500 Lithium-ion Battery 3000 mAh



### Product overview: Fluke BP500 Lithium-ion Battery 3000 mAh

Replacement or spare high-capacity 7.4 V 3000MAH Li-ion battery for Fluke 500 Series Battery Analyzers.

## Ordering information



### **Fluke BP500**

Lithium-ion Battery 3000 mAh

---

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**Fluke Australia**  
Unit 26, 7 Anella Ave  
Castle Hill, NSW 2154 Australia  
Phone: 61 2 8850-3333  
[www.fluke.com.au](http://www.fluke.com.au)

©2021 Fluke Corporation. All rights reserved.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**