

TECHNICAL DATA

Fluke BTL20 Interactive Battery Analyzer Test Probe



Product overview: Fluke BTL20 Interactive Battery Analyzer Test Probe

Intelligent test probe set for Fluke BT520 Battery Analyzer (temperature measurement function not included). The BTL20 has a built in LCD display, tri-light LED for pass/fail warning indication, voice count for battery ID and integrated save button. Comes with both long and short probe extenders for use in hard to reach battery cabinets.

Ordering information



Fluke BTL20

Interactive Battery Analyzer Test Probe

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

Fluke Australia
Unit 26, 7 Anella Ave
Castle Hill, NSW 2154 Australia
Phone: 61 2 8850-3333
www.fluke.com.au

©2021 Fluke Corporation. All rights reserved.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**