

TECHNICAL DATA

Fluke RS200 Probe-Accessories Replacement Set



Key features

Compatible with PM8918, VPS100, VPS101, VPS121 and VPS200-series probes.

Product overview: Fluke RS200 Probe-Accessories Replacement Set

Contains replacement parts for usage with the VPS200 & PM8918 Scope Probes.

Ordering information



Fluke RS200

ul>

- Set includes 600V CAT III / 1000V CAT II rated:
 - Hook Clip for Probe Tip (2 Red, 1 Gray)
 - 4-mm Test Probe for Probe (1 Red, 1 Gray)
 - Ground Lead with Hook Clip (1 Red, 1 Gray)
 - Ground Lead with Mini Alligator (2 Black).
 - Ground Spring for Probe Tip (5x) 30V CAT II Rated
- Max. Input Voltage: 600V CAT III, 1000V CAT II per IEC 1010-2-031
- Max. Voltage to Ground: 600V CAT III, 1000V CAT II per IEC 1010-2-031

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

Fluke Australia
Unit 26, 7 Anella Ave
Castle Hill, NSW 2154 Australia
Phone: 61 2 8850-3333
www.fluke.com.au

©2021 Fluke Corporation. All rights reserved.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**