

Fluke TLK290 Flexible Socket Probe Set



Key features

- Test probe set to be used on motor and three phase sockets
- Probes have flexible test points that fits securely between 4 and 8 mm sockets
- For use with 4 mm shrouded test leads, such as Fluke TL224
- Safe contact e.g. in CEE 16 A und CEE 32 A plugs
- Probe CAT II 1000V, 8A; Alligator CAT III 1000V, 16A

Product overview: Fluke TLK290 Flexible Socket Probe Set

Set includes three Flexible Socket Probes and a Large Alligator Clip.

Ordering information



Fluke TLK290



Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

Fluke Australia
Unit 26, 7 Anella Ave
Castle Hill, NSW 2154 Australia
Phone: 61 2 8850-3333
www.fluke.com.au

©2021 Fluke Corporation. All rights reserved.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**