

# Fluke TL221 SureGrip™ Test Lead Extenders



## Key features

- Extends test leads by 1.5 meters
- Silicone-insulated wire resists heat & cold
- Accepts wide range of Fluke test probes and clips

## Product overview: Fluke TL221 SureGrip™ Test Lead Extenders

The Fluke TL221 SureGrip™ Test Lead Extenders includes silicone leads and lead couplers with reinforced strain relief. With 1.5 meters of added length you can reach test points that were previously inaccessible. Rubber over-molded surfaces and finger-hugging curves gives you a comfortable, reliable grip.

### Designed for safety

The Fluke TL221 SureGrip™ Test Lead Extenders are CAT III 1000 V, CAT IV 600 V rated.

### Compatible with Fluke test probes and clips

TL221 lead extensions are compatible with a wide range of Fluke test probes and clips. Some often-used probe tips include:

- TP175 TwistGuard™ Test Probes

- TP220 SureGrip™ Industrial Test Probes
- AC220 SureGrip™ Alligator Clips
- AC285 SureGrip™ Alligator Clips
- AC280 SureGrip™ Hook Clips
- TP1 SlimReach™ Test Probes
- AC283 SureGrip™ Pincer Clips

## Ordering information



### **Fluke TL221**

- One pair (red, black) of silicone insulated leads with straight connectors on both ends
- Reinforced strain relief
- Includes one pair (red, black) of female couplers
- Extends test leads by 1.5 m
- 600 V CAT IV, 1000 V CAT III, 10 A rating. UL listed
- One year warranty

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**Fluke Australia**  
Unit 26, 7 Anella Ave  
Castle Hill, NSW 2154 Australia  
Phone: 61 2 8850-3333  
[www.fluke.com.au](http://www.fluke.com.au)

©2021 Fluke Corporation. All rights reserved.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**