

# Fluke TLK287 Electronics Master Test Lead Set



## Key features

Precision electronic probes with multiple spring-loaded sharp tips to maximize contact with SMD test points. Micro-grabbers and leads to enable fine SMD test. Mini-alligator clips, grabbers, and probes for virtually every electronics need.

## Product overview: Fluke TLK287 Electronics Master Test Lead Set

Fluke TLK287 Electronics Master Test Lead Set is ideal for electronic test and design. The set includes precision electronic probes with multiple spring-loaded sharp tips to maximize contact with SMD test points. Micro-grabbers and leads enable fine SMD testing. Mini-alligator clips, grabbers and probes for virtually every electronics need. And modular test leads and lead-probe combinations for versatility and couplers for long lead applications.

## Ordering information



### **Fluke-TLK287**

---

Includes:

- Precision electronic probes with multiple spring-loaded sharp tips to maximize contact with SMD test points
  - Micro-grabbers and leads to enable fine SMD test
  - Mini-alligator clips, grabbers, and probes for virtually every electronics need
-



**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**Fluke Australia**  
Unit 26, 7 Anella Ave  
Castle Hill, NSW 2154 Australia  
Phone: 61 2 8850-3333  
[www.fluke.com.au](http://www.fluke.com.au)

©2021 Fluke Corporation. All rights reserved.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**