

Fluke Networks MT-8203-20 IntelliTone Test Leads



Key features

- Test leads with Bed of Nails for Pro 100/200
- Two Banana Jack test leads with piercing pin and alligator clips
- U-shaped nose for blade-shaped terminals
- Notched jaw and serrated teeth for gripping
- Large piercing needle for large gauge wires

Product overview: Fluke Networks MT-8203-20 IntelliTone Test Leads

IntelliTone™ MT-8203-20 test leads featuring an insulation-piercing clip on one end and a shrouded, 4-mm right-angle banana plug on the other. The clip features a bed-of-nails for piercing 22-28 AWG wires.

Ordering information



MT-8203-20

Fluke Networks MT-8203-20 IntelliTone™ Test Leads

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

Fluke Australia
Unit 26, 7 Anella Ave
Castle Hill, NSW 2154 Australia
Phone: 61 2 8850-3333
www.fluke.com.au

©2021 Fluke Corporation. All rights reserved.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**