

TECHNICAL DATA

Fluke 700M20TH Premium Transmitter Test Hose Kit



Product overview: Fluke 700M20TH Premium Transmitter Test Hose Kit

Easy to use premium pressure transmitter test hose kit for use with female M20 connections

- 5,000 psi rating
- Includes 1/8" NPT connection for permanent mount to a calibrator or test pump
- No tools required, male connector for connection to female M20 connections

Ordering information



Fluke 700M20TH

Includes:

- One (1) 5,000 psi rated, premium test hose
 - Premium, no tools required M20 adapter
-

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com/en-gb

©2022 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
01/2022

**Modification of this document is not permitted
without written permission from Fluke Corporation.**