

## TECHNICAL DATA

# Fluke 700MTH Premium Transmitter Test Hose Kit



## Product overview: Fluke 700MTH Premium Transmitter Test Hose Kit

Easy to use premium pressure transmitter test hose kit for connection to 1/4" female metric/BSP

- 5,000 psi rating
- Includes 1/8" NPT connection for permanent mount to a calibrator or test pump
- No tools required, male connector for connection to 1/4" female metric/BSP

## Ordering information



### Fluke 700MTH

---

Includes:

- One (1) 5,000 psi rated, premium test hose
  - Premium, no tools required male ¼" metric/BSP connector
-

**Fluke.** *Keeping your world up and running.®*

**Fluke (UK) Ltd.**  
52 Hurricane Way  
Norwich, Norfolk  
NR6 6JB  
United Kingdom  
Tel.: +44 (0)20 7942 0708  
E-mail: [cs.uk@fluke.com](mailto:cs.uk@fluke.com)  
[www.fluke.com/en-gb](http://www.fluke.com/en-gb)

©2022 Fluke Corporation. All rights reserved.  
Data subject to alteration without notice.  
01/2022

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**