


TECHNICAL DATA

Fluke ES-162P4-2 4-Pole Stake Kit



Product overview: Fluke ES-162P4-2 4-Pole Stake Kit

Needed to perform 3- and 4-Pole Fall-of-Potential testing (Stake testing). Used with both the Fluke 1623-2 and Fluke 1625-2. Already included in the Fluke 1625-2 Kit.

Ordering information



Fluke ES-162P4-2

- (4) Stakes
- (1) 50 m cable reel of red wire
- (1) 25 m cable reel of blue wire
- (1) 25 m cable reel of green wire

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com/en-gb

©2022 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
01/2022

**Modification of this document is not permitted
without written permission from Fluke Corporation.**