

TECHNICAL DATA

TLK281 -1, Automotive Test Lead Kit

Product overview: TLK281 -1, Automotive Test Lead Kit

Complete set of test and measurement accessories designed especially for automotive specialists. Use the sharp stainless steel probes to contact the wire conductor or pinpoint tests on a circuit or relay, and the large jaw alligator clips to get a good ground on the engine or to connect to the battery terminals. An extra long set of alligator clips easily grip to the alternator BAT terminal or the S terminals. The set comes with silicone test leads and a case to keep it organized.

- Tested to 60 V DC

Kit contains:

- **TP81 Insulation Piercing Clip Set**
- **TL224 Suregrip™ Silicone Test Lead Set**
- **TP220-1 Suregrip™ Test Probe Set**
- **AC220 Suregrip™ Alligator Clip Set**
- **AC285 Suregrip™ Large Jaw Alligator Clip Set**
- **Soft Case**

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com/en-gb

©2022 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
01/2022

**Modification of this document is not permitted
without written permission from Fluke Corporation.**