

TECHNICAL DATA


Fluke Networks Replacement Remote ID #1 / Wiremapper for LinkIQ™



Product overview: Fluke Networks Replacement Remote ID #1 / Wiremapper for LinkIQ™

Replacement Remote ID #1 / Wiremapper for [LinkIQ™ Cable+Network Tester](#). Works with the [Industrial Ethernet version](#) of LinkIQ as well.

Ordering information



REMOTEID-1

Replacement Remote ID #1 / Wiremapper for LinkIQ™

Fluke. Keeping your world up and running.®

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

BUT. FLUKE SOUTH EAST ASIA PTE LTD
Menera Satu Sentra Kelapa Gading #06-05
Jl. Bulevar Kelapa Gading Kav. LA# No. 1
Summarecon Kelapa Gading
Jakarta Utara 14240
Indonesia
Tel: +62 21 2938 5922
Email: info.asean@fluke.com
www.fluke.com/id-en

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500

©2021 Fluke Corporation. Specifications subject to
change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**