

TECHNICAL DATA

Fluke 80PK-1 Bead Probe



Key features

- Measurement range: -40 to 260 °C
- One year warranty

Product overview: Fluke 80PK-1 Bead Probe

Type-K thermocouple for general purpose applications

Ordering information



Fluke 80PK 1

PROBE, THERMOCOUPLE, BEADED K-TYPE

- Probe Length: 0.04 IN
 - Lead Length: 39 IN (25.40 MM)
 - Product Height: 1.50 IN
-

Fluke. *Keeping your world up and running.*®

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa

+31 (0)40 267 5100

In Canada (800)-36-FLUKE

From other countries +1 (425) 446-5500

BUT. FLUKE SOUTH EAST ASIA PTE LTD

Menera Satu Sentra Kelapa Gading #06-05

Jl. Bulevar Kelapa Gading Kav. LA# No. 1

Summarecon Kelapa Gading

Jakarta Utara 14240

Indonesia

Tel: +62 21 2938 5922

Email: info.asean@fluke.com

www.fluke.com/id-en

©2022 Fluke Corporation. Specifications subject to change without notice.

01/2022

Modification of this document is not permitted without written permission from Fluke Corporation.