

Fluke TLK290 Flexible Socket Probe Set



Key features

- Test probe set to be used on motor and three phase sockets
- Probes have flexible test points that fits securely between 4 and 8 mm sockets
- For use with 4 mm shrouded test leads, such as Fluke TL224
- Safe contact e.g. in CEE 16 A und CEE 32 A plugs
- Probe CAT II 1000V, 8A; Alligator CAT III 1000V, 16A

Product overview: Fluke TLK290 Flexible Socket Probe Set

Set includes three Flexible Socket Probes and a Large Alligator Clip.

Ordering information



Fluke TLK290

Fluke. *Keeping your world up and running.*®

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500

BUT. FLUKE SOUTH EAST ASIA PTE LTD
Menera Satu Sentra Kelapa Gading #06-05
JI. Bulevar Kelapa Gading Kav. LA# No. 1
Summarecon Kelapa Gading
Jakarta Utara 14240
Indonesia
Tel: +62 21 2938 5922
Email: info.asean@fluke.com
www.fluke.com/id-en
©2021 Fluke Corporation. Specifications subject to
change without notice.
10/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**