

Fluke TLK-220 US SureGrip™ Industrial Test Lead



Key features

Handy all-in-one kit that includes AC220 SureGrip Alligator Clip Set, AC285 SureGrip Large Jaw Alligator Clip Set, TP220 Sharp Test Probes, TL224 SureGrip Silicone Light to Straight Test Leads, Zippered Vinyl Carry Case with Moveable Divider.

Product overview: Fluke TLK-220 US SureGrip™ Industrial Test Lead

Special value SureGrip Set for industrial applications. Includes two sizes of alligator clips and sharp test probes for reliable connections, test leads and a roomy carry case to hold and protect your meter and all your accessories.

Ordering information



Fluke TLK 220

US SureGrip™ Industrial Test Lead

Included with set:

- AC220 Plunger Style Alligator Clips
 - AC285 Large Jaw Alligator Clips
 - TP220 Sharp Test Probes
 - TL224 Right to Straight Test Leads
 - Zippered Vinyl Carry Case with Moveable Divider
-

Fluke. *Keeping your world up and running.*®

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500

BUT. FLUKE SOUTH EAST ASIA PTE LTD
Menera Satu Sentra Kelapa Gading #06-05
JI. Bulevar Kelapa Gading Kav. LA# No. 1
Summarecon Kelapa Gading
Jakarta Utara 14240
Indonesia
Tel: +62 21 2938 5922
Email: info.asean@fluke.com
www.fluke.com/id-en
©2021 Fluke Corporation. Specifications subject to
change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**