

TECHNICAL DATA

Fluke C43 Soft Carrying Case



Key features

- Zippered carrying case with durable polyester exterior
- Inner front pocket with closure and two pockets to store test leads or small accessories
- Includes 1 detachable 10 inch handle and 1 detachable 42 inch shoulder strap
- Companion case to Fluke 1520 Tester and One Touch Network Analyzers
- Also ideal for Fluke 40 series power meters and i1000s, i3000s large current clamps
- External dimensions: 32 x 19 x 9 cm (12.5 x 7.5 x 3.5 in)

Ordering information



Fluke C43

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com

©2021 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**