

TECHNICAL DATA

Fluke 700TTH10K Premium Transmitter Test Hose Kit



Product overview: Fluke 700TTH10K Premium Transmitter Test Hose Kit

Easy to use premium pressure transmitter test hose kits for use with a calibrator or test pump.

- 10,000 psi rating
- Included 1/8" NPT connection for permanent mount to a calibrator or test pump
- No tools required, 1/4" NPT male connector for connection to a transmitter

Ordering information



Fluke 700TTH10K

Includes:

- One (1) premium hose kit and adapters
 - Premium/no tools required 1/4" NPT transmitter connection fitting
-

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com

©2021 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**