

TECHNICAL DATA

Fluke 80PK-1 Bead Probe



Key features

- Measurement range: -40 to 260 °C
- One year warranty

Product overview: Fluke 80PK-1 Bead Probe

Type-K thermocouple for general purpose applications

Ordering information



Fluke 80PK 1

PROBE, THERMOCOUPLE, BEADED K-TYPE

- Probe Length: 0.04 IN
 - Lead Length: 39 IN (25.40 MM)
 - Product Height: 1.50 IN
-

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com

©2021 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**