

TECHNICAL DATA

Fluke 80PK-25 SureGrip™ Piercing Temperature Probe



Key features

- Type-K thermocouple suitable for food industry, liquids and gels
- Measurement range : -40 to 260 °C
- One year warranty

Ordering information



Fluke 80PK 25

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com

©2021 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**