


# Fluke AS200-R Probe Accessory Set, Red




## Product overview: Fluke AS200-R Probe Accessory Set, Red

Probe Accessory Set for usage with VPS200 & PM8918 Series Scope Probes.

## Specifications: Fluke AS200-R Probe Accessory Set, Red

Specifications	
Max. input voltage	600 V CAT III, 1000V CAT II per IEC 1010-2-031
Max. voltage to ground	600 V CAT III, 1000V CAT II per IEC 1010-2-031
Compatible with	VPS100, VPS101, VPS121, VPS200-Series and PM8918-Series probes

## Ordering information



### **Fluke AS200 R**

Set includes 600V CAT III / 1000V CAT II rated:

- Industrial Alligator for Probe
- 2-mm Test Probe for Probe
- 2-mm Test Probe for Banana
- Industrial Alligator for Banana
- Ground Lead 4-mm Banana

**Fluke.** *Keeping your world up and running.®*

**Fluke (UK) Ltd.**  
52 Hurricane Way  
Norwich, Norfolk  
NR6 6JB  
United Kingdom  
Tel.: +44 (0)20 7942 0708  
E-mail: [cs.uk@fluke.com](mailto:cs.uk@fluke.com)  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation. All rights reserved.  
Data subject to alteration without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**