

TECHNICAL DATA

Fluke RS200 Probe-Accessories Replacement Set




Key features

Compatible with PM8918, VPS100, VPS101, VPS121 and VPS200-series probes.

Product overview: Fluke RS200 Probe-Accessories Replacement Set

Contains replacement parts for usage with the VPS200 & PM8918 Scope Probes.

Ordering information



Fluke RS200

ul>

- Set includes 600V CAT III / 1000V CAT II rated:
 - Hook Clip for Probe Tip (2 Red, 1 Gray)
 - 4-mm Test Probe for Probe (1 Red, 1 Gray)
 - Ground Lead with Hook Clip (1 Red, 1 Gray)
 - Ground Lead with Mini Alligator (2 Black).
 - Ground Spring for Probe Tip (5x) 30V CAT II Rated
- Max. Input Voltage: 600V CAT III, 1000V CAT II per IEC 1010-2-031
- Max. Voltage to Ground: 600V CAT III, 1000V CAT II per IEC 1010-2-031

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com

©2021 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**