

TECHNICAL DATA

Fluke TL970 Hook and Pincer Set



Key features

Handy all-in-one kit that includes 3 Test Lead Sets:

- TL940 Mini-Hook Test Lead Set
- TL950 Mini-Pincer Test Lead Set
- TL960 Micro-Hook Test Lead Set
- One year warranty

Ordering information



Fluke TL970

Handy all-in-one kit that includes:

- 3 Test Lead Sets:
 - TL940 Mini-Hook Test Lead Set
 - TL950 Mini-Pincer Test Lead Set
 - TL960 Micro-Hook Test Lead Set
-

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com

©2021 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**