

Fluke TLK-225 SureGrip™ Master Accessory Set



Key features

The Fluke TLK-225 SureGrip™ Master Accessory Set is the perfect replacement kit that provides flexibility and comfort with all SureGrip™ leads and probes.

Set includes:

- TL224 Test Lead Set
- TP175 TwistGuard™ Test Probes
- AC285 SureGrip™ Alligator Clips
- AC220 SureGrip™ Alligator Clips
- AC280 SureGrip™ Hook Clips
- AC283 SureGrip™ Pincer Clips
- A convenient six-pocket roll up pouch
- CAT III 1000 V, CAT IV 600 V safety rating

Product overview: Fluke TLK-225 SureGrip™ Master Accessory Set

Test leads are an integral part of the complete measurement system. Worn out test leads can cause inaccurate readings and pose serious shock or electrocution hazard. The Fluke TLK-225 SureGrip™ Master Accessory Set with 13 items is the perfect replacement kit that provides flexibility and comfort with all SureGrip™ leads and probes in a convenient six pocket, roll up pouch.

Ordering information



Fluke TLK225

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com

©2021 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**