

## Fluke VPS101 – Voltage Probe Set, 1:1, one black



### Key features

- 300 V CAT III rated
  - 1:1 Voltage Probe
  - Hook Clip for Probe Tip
  - Ground Lead with Mini Alligator
  - Insulation Sleeve
- 30 V CAT II rated Ground Spring for Probe Tip

### Product overview: Fluke VPS101 – Voltage Probe Set, 1:1, one black

Double insulated voltage probe for ScopeMeter 190 series.

### Specifications: Fluke VPS101 – Voltage Probe Set, 1:1, one black

Specification	
Attenuation	1x
Bandwidth	30 MHz
Max. input voltage	300 V CAT III
Max. voltage to ground	600 V CAT III per IEC / EN61010-031



Input resistance	1 M $\Omega$
Input capacitance	42 pF + input capacitance of instrument
Cable length	1.2 m
Warranty	One year

**Fluke.** *Keeping your world up and running.®*

**Fluke (UK) Ltd.**  
52 Hurricane Way  
Norwich, Norfolk  
NR6 6JB  
United Kingdom  
Tel.: +44 (0)20 7942 0708  
E-mail: [cs.uk@fluke.com](mailto:cs.uk@fluke.com)  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation. All rights reserved.  
Data subject to alteration without notice.  
11/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**