

VPS220-G - Voltage Probe Set, 100:1, one grey



Key features

Each probe set includes:

- 600 V CAT III / 1000 V CAT II rated 100:1 Voltage Probe
- Hook Clip for Probe Tip
- Ground Lead with Mini Alligator
- Insulation Sleeve
- 30 V CAT II rated Ground Spring for Probe Tip

Product overview: VPS220-G - Voltage Probe Set, 100:1, one grey


Double insulated voltage probe for ScopeMeter 190 series.

Specifications: VPS220-G - Voltage Probe Set, 100:1, one grey

Specification	
Attenuation	100:1
Bandwidth	200 MHz
Max. input voltage	1000 V CAT II or 600 V CAT III

Max. voltage to ground	1000 V CAT II or 600 V CAT III per IEC / EN61010-031
Input resistance	100 MΩ
Input capacitance	4.6 pF
Compensation range	10 to 25 pF
Cable length	1.2 m

Ordering information



Fluke VPS220-G

Voltage Probe Set, 100:1, one grey

Each probe set includes:

- 600 V CAT III / 1000 V CAT II rated 100:1 Voltage Probe
 - Hook Clip for Probe Tip
 - Ground Lead with Mini Alligator
 - Insulation Sleeve
 - 30 V CAT II rated Ground Spring for Probe Tip
-

Fluke. *Keeping your world up and running.®*

Fluke (UK) Ltd.
52 Hurricane Way
Norwich, Norfolk
NR6 6JB
United Kingdom
Tel.: +44 (0)20 7942 0708
E-mail: cs.uk@fluke.com
www.fluke.com

©2021 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**