

# Fluke Calibration P5521 and P5522 Liquid-to-Liquid Separators



## Key features

- P5521-For pressures to 10,000 psi (700 bar)
- P5522-For pressures to 7,000 psi (480 bar) for Skydrol or other aggressive fluids

## Product overview: Fluke Calibration P5521 and P5522 Liquid-to-Liquid Separators

These liquid-to-liquid separators connect directly to the test port of a hydraulic deadweight tester or comparison test pump. A flexible diaphragm separates the fluids, protecting the calibrator from contamination, and allows calibration of the device in its specific working fluid.

## Specifications: Fluke Calibration P5521 and P5522 Liquid-to-Liquid Separators

Specifications	
Pressure Range	5521: 0 to 10,000 psi (700 bar) 5522: 0 to 7,000 psi (480 bar) for Skydrol and other aggressive fluids
Instrument Size (W x H)	4 x 6 in (100 x 150 mm)



Instrument Weight	7.5 lb (3.5 kg)
Test Port Adapter	Uses adapters supplied with calibrator
Calibrator Connection	3/8 BSP male (mating adapter supplied)
Body Material	Aluminum bronze
Diaphragm	5521- Viton® 5522-Ethylene Propylene
Reservoir Material	Hard-anodized aluminum

## Ordering information



### **P5521 and P5522**

5521 and 5522 Liquid to Liquid Separators

---

Liquid-to-Liquid Separator

Pressure Range 0 to 10,000 psi (700 bar)

---

**Fluke.** *Keeping your world up and running.®*

**Fluke (UK) Ltd.**  
52 Hurricane Way  
Norwich, Norfolk  
NR6 6JB  
United Kingdom  
Tel.: +44 (0)20 7942 0708  
E-mail: [cs.uk@fluke.com](mailto:cs.uk@fluke.com)  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation. All rights reserved.  
Data subject to alteration without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**