

TECHNICAL DATA

## Fluke H6 Infrared Thermometer Holster



### Key features

- Durable nylon construction protects your thermometer while keeping it within reach
- Cuff hook attachment closure provide quick access and a secure closure
- Extra accessory pocket for storing temperature probe
- Snap-open loop attaches easily to your belt

### Product overview: Fluke H6 Infrared Thermometer Holster

Durable nylon construction protects your thermometer while keeping it within reach, cuff hook attachment closure provide quick access and a secure closure.

## Ordering information



### **Fluke H6**

Fluke H6 Infrared Thermometer Holster

---



**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Canada (800) 36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**