

TECHNICAL DATA

Fluke 80PK-22 SureGrip™ Immersion Temperature Probe



Key features

- Type-K thermocouple for use in liquids and gels
- Measurement range: -40 to 1090 °C
- One year warranty

Product overview: Fluke 80PK-22 SureGrip™ Immersion Temperature Probe

Fluke 80PK-22 SureGrip™ Immersion Temperature Probe is designed for use in liquid and gels. It can also be used as a general purpose temperature probe. The thermocouple junction is protected from tip to handle by an Inconel sheath. The 80PK-22 can be used with any temperature-measuring instrument that is designed to accept type K thermocouples and has a miniature connector input.

Ordering information



Fluke 80PK 22

Fluke 80PK-22 SureGrip™ Immersion Temperature Probe

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**