

TECHNICAL DATA

Fluke TL910 Electronic Test Probes



Key features

- Very small 1 mm tips to access hard to reach test points
- Probe tip length up to 100 mm, test lead length: 1 m
- 3 sets of spring loaded gold tips and 2 sets of stainless steel tips included
- CAT II 1000 V, 3 A rating
- One year warranty

Ordering information



Fluke TL910



Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**