

TECHNICAL DATA

## Fluke C43 Soft Carrying Case



### Key features

- Zippered carrying case with durable polyester exterior
- Inner front pocket with closure and two pockets to store test leads or small accessories
- Includes 1 detachable 10 inch handle and 1 detachable 42 inch shoulder strap
- Companion case to Fluke 1520 Tester and One Touch Network Analyzers
- Also ideal for Fluke 40 series power meters and i1000s, i3000s large current clamps
- External dimensions: 32 x 19 x 9 cm (12.5 x 7.5 x 3.5 in)

## Ordering information



**Fluke C43**

---

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**

PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa +31 (0)40 267 5100

In Canada (800)-36-FLUKE

From other countries +1 (425) 446-5500

**Fluke Electronics (Malaysia) Sdn Bhd**

Suite D-3A-4, Block D

SouthGage Commercial Centre, No2

Jalan Dua Off Jln Chan Sow Lin

55200 KL

Malaysia

Email: [info.asean@fluke.com](mailto:info.asean@fluke.com)

[www.fluke.com/my](http://www.fluke.com/my)

©2021 Fluke Corporation. Specifications subject to change without notice.

12/2021

**Modification of this document is not permitted without written permission from Fluke Corporation.**