

Fluke 700PTPK Pneumatic Test Pressure Kit



Key features

- Easy to use, rugged and reliable construction, capable of generating pneumatic pressure from -12.7 psi/ -0.87 bar to 600 psi/41 bar
- Combine with one of the Fluke 700G Precision Pressure Gauges for a complete pressure testing and calibration solution
- In addition, combine with the Fluke 700G/TRACK Software to upload over 8,000 pressure measurements logged in the field to a PC

Product overview: Fluke 700PTPK Pneumatic Test Pressure Kit

Generates pressure up to 600 psi, 41 bar. Includes the 700PTP-1 pneumatic hand pump, test hoses and adapters, and protective hard case.

Fluke. *Keeping your world up and running.®*

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa +31 (0)40 267 5100

In Canada (800)-36-FLUKE

From other countries +1 (425) 446-5500

Fluke Electronics (Malaysia) Sdn Bhd

Suite D-3A-4, Block D

SouthGage Commercial Centre, No2

Jalan Dua Off Jln Chan Sow Lin

55200 KL

Malaysia

Email: info.asean@fluke.com

www.fluke.com/my

©2021 Fluke Corporation. Specifications subject to change without notice.

11/2021

Modification of this document is not permitted without written permission from Fluke Corporation.