


TECHNICAL DATA

Fluke ES-162P4-2 4-Pole Stake Kit



Product overview: Fluke ES-162P4-2 4-Pole Stake Kit

Needed to perform 3- and 4-Pole Fall-of-Potential testing (Stake testing). Used with both the Fluke 1623-2 and Fluke 1625-2. Already included in the Fluke 1625-2 Kit.

Ordering information



Fluke ES-162P4-2

- (4) Stakes
- (1) 50 m cable reel of red wire
- (1) 25 m cable reel of blue wire
- (1) 25 m cable reel of green wire

Fluke. *Keeping your world up and running.®*

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa +31 (0)40 267 5100

In Canada (800)-36-FLUKE

From other countries +1 (425) 446-5500

Fluke Electronics (Malaysia) Sdn Bhd

Suite D-3A-4, Block D

SouthGage Commercial Centre, No2

Jalan Dua Off Jln Chan Sow Lin

55200 KL

Malaysia

Email: info.asean@fluke.com

www.fluke.com/my

©2021 Fluke Corporation. Specifications subject to change without notice.

11/2021

Modification of this document is not permitted without written permission from Fluke Corporation.