

# Fluke TLK-220 US SureGrip™ Industrial Test Lead



## Key features

Handy all-in-one kit that includes AC220 SureGrip Alligator Clip Set, AC285 SureGrip Large Jaw Alligator Clip Set, TP220 Sharp Test Probes, TL224 SureGrip Silicone Light to Straight Test Leads, Zippered Vinyl Carry Case with Moveable Divider.

## Product overview: Fluke TLK-220 US SureGrip™ Industrial Test Lead

Special value SureGrip Set for industrial applications. Includes two sizes of alligator clips and sharp test probes for reliable connections, test leads and a roomy carry case to hold and protect your meter and all your accessories.

## Ordering information



### **Fluke TLK 220**

US SureGrip™ Industrial Test Lead

---

Included with set:

- AC220 Plunger Style Alligator Clips
  - AC285 Large Jaw Alligator Clips
  - TP220 Sharp Test Probes
  - TL224 Right to Straight Test Leads
  - Zippered Vinyl Carry Case with Moveable Divider
-

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Europe/M-East/Africa +31 (0)40 267 5100  
In Canada (800)-36-FLUKE  
From other countries +1 (425) 446-5500

**Fluke Electronics (Malaysia) Sdn Bhd**  
Suite D-3A-4, Block D  
SouthGage Commercial Centre, No2  
Jalan Dua Off Jln Chan Sow Lin  
55200 KL  
Malaysia  
Email: [info.asean@fluke.com](mailto:info.asean@fluke.com)  
[www.fluke.com/my](http://www.fluke.com/my)

©2021 Fluke Corporation. Specifications subject to  
change without notice.  
11/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**