

Fluke VPS101 – Voltage Probe Set, 1:1, one black



Key features

- 300 V CAT III rated
 - 1:1 Voltage Probe
 - Hook Clip for Probe Tip
 - Ground Lead with Mini Alligator
 - Insulation Sleeve
- 30 V CAT II rated Ground Spring for Probe Tip

Product overview: Fluke VPS101 – Voltage Probe Set, 1:1, one black

Double insulated voltage probe for ScopeMeter 190 series.

Specifications: Fluke VPS101 – Voltage Probe Set, 1:1, one black

Specification	
Attenuation	1x
Bandwidth	30 MHz
Max. input voltage	300 V CAT III
Max. voltage to ground	600 V CAT III per IEC / EN61010-031

Input resistance	1 MΩ
Input capacitance	42 pF + input capacitance of instrument
Cable length	1.2 m
Warranty	One year

Fluke. *Keeping your world up and running.®*

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa +31 (0)40 267 5100

In Canada (800)-36-FLUKE

From other countries +1 (425) 446-5500

Fluke Electronics (Malaysia) Sdn Bhd

Suite D-3A-4, Block D

SouthGage Commercial Centre, No2

Jalan Dua Off Jln Chan Sow Lin

55200 KL

Malaysia

Email: info.asean@fluke.com

www.fluke.com/my

©2021 Fluke Corporation. Specifications subject to change without notice.

12/2021

Modification of this document is not permitted without written permission from Fluke Corporation.