

TECHNICAL DATA


Fluke Networks Replacement Remote ID #1 / Wiremapper for LinkIQ™



Product overview: Fluke Networks Replacement Remote ID #1 / Wiremapper for LinkIQ™

Replacement Remote ID #1 / Wiremapper for [LinkIQ™ Cable+Network Tester](#). Works with the [Industrial Ethernet version](#) of LinkIQ as well.

Ordering information



REMOTED-1

Replacement Remote ID #1 / Wiremapper for LinkIQ™

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**