

TECHNICAL DATA

Fluke 700MTH Premium Transmitter Test Hose Kit



Product overview: Fluke 700MTH Premium Transmitter Test Hose Kit

Easy to use premium pressure transmitter test hose kit for connection to 1/4" female metric/BSP

- 5,000 psi rating
- Includes 1/8" NPT connection for permanent mount to a calibrator or test pump
- No tools required, male connector for connection to 1/4" female metric/BSP

Ordering information



Fluke 700MTH

Includes:

- One (1) 5,000 psi rated, premium test hose
 - Premium, no tools required male 1/4" metric/BSP connector
-

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**