

TECHNICAL DATA

# Fluke L210 Probe Light Kit




## Key features

- Bright white LED illuminates contact area
- 8" probe extenders keep hands away from live circuits
- Extender fits between modular test probe and test lead
- Small, rugged light frees both hands for work
- 120 hours of battery life
- Uses 2 3-volt watch batteries

## Product overview: Fluke L210 Probe Light Kit

L210 Probe Light and Probe Extenders

## Ordering information



### **Fluke L210**

Includes L200 Probe Light and TP280 Test Probe Extenders



**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Canada (800) 36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**