


TECHNICAL DATA

Fluke PM9094/101 Mini Testhook Set for PM8918 series probes




Product overview: Fluke PM9094/101 Mini Testhook Set for PM8918 series probes

Set consists of: 4 Mini Testhooks, 2 red, 2 gray. Max. voltage: 30 V rms.

Note: When using the mini Test-hooks any potential must not exceed 42 Vpk/30 Vrms

Ordering information



Fluke PM9094/101

Mini Testhook Set for PM8918 series probes

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**