

TECHNICAL DATA

Fluke 700 PMP Pressure Pump



Product overview: Fluke 700 PMP Pressure Pump

Description: The Fluke 700PMP is a hand-operated pressure pump to provide pressures up to 150 psi/1000 kPa. Output fitting is 1/8 FNPT.

Application: Linear stroke of 4 cm. Multi-turn vernier for fine adjustment of pressure.

Specifications: Fluke 700 PMP Pressure Pump

Mechanical and General Specifications	
Size	265 x 35 mm
Weight	340 g
Warranty	One year

Ordering information



Fluke 700PMP

PMP Pressure Pump

Instruction sheet



Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**