

TECHNICAL DATA

Fluke 376 Clamp / PRV240 Proving Unit Kit



Ordering information



FLUKE-376-PRV240

Fluke 376 Clamp / PRV240 Proving Unit Kit

Includes:

True-rms AC/DC Clamp Meter with iFlex™

- 18-inch iFlex™ flexible current probe
- Coated instruction card
- Safety information sheet
- Soft carrying case
- TL75 test leads
- Two AA Alkaline batteries (installed)

Fluke PRV240 Proving Unit

- TPAK hanging strap
- Four AA batteries
- Instruction sheet



Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2022 Fluke Corporation.
Specifications subject to change without notice.
01/2022

**Modification of this document is not permitted
without written permission from Fluke Corporation.**