

TECHNICAL DATA

Fluke H6 Infrared Thermometer Holster



Key features

- Durable nylon construction protects your thermometer while keeping it within reach
- Cuff hook attachment closure provide quick access and a secure closure
- Extra accessory pocket for storing temperature probe
- Snap-open loop attaches easily to your belt

Product overview: Fluke H6 Infrared Thermometer Holster

Durable nylon construction protects your thermometer while keeping it within reach, cuff hook attachment closure provide quick access and a secure closure.

Ordering information



Fluke H6

Fluke H6 Infrared Thermometer Holster

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**