

TECHNICAL DATA

Fluke 80PJ-1 Bead Probe



Product overview: Fluke 80PJ-1 Bead Probe

- Type-J thermocouple for general purpose applications
- Measurement range: -40 to 260°C
- One year warranty

Ordering information



Fluke 80PJ-1

Bead Probe

-
- Type-K thermocouple for general purpose applications
 - Measurement range: -40 to 260 °C
 - One year warranty
-



Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**